

SuperBall®

Meniscal Repair System



Complete Solution for All-inside Meniscal Repairs

SuperBall, the all-inside, all-suture, knotless meniscal repair system, now provides you with a complete solution to access and repair even the hardest to reach areas of the meniscus. Our two newest curvatures utilize the same great technology you have come to trust in maintaining meniscal form and function.

- **Strong and reliable** repair using all-suture implants
- **Proprietary suture mesh** that lays flat and reduces tissue cut-through
- Unique integrated **tensioning mechanism** that mimics inside-out tensioning and reduces anchor pull-outs




12° Reverse




24° Up

Multiple Delivery Options



18° Up: Standard design for majority of tears



24° Up: Designed for hard-to-reach posterior tears and more anterior midbody tears



12° Reverse: Designed for under-surface tears



ArcuroMedical.com